



DAC-SPECTROMED S.R.L.

МД-2025, Молдова, г. Кишинев, Ул. Гестемицану, 20
Тел.: /+3732/ 739961, 739968; факс: /+3732/ 727522
Email: office@dacspectromed.com

DIAGNOSTIC REAGENT KITS FOR CLINICAL BIOCHEMISTRY

Product/Description/Method	Total reagent Vol., mL
ACP - DAC ACID PHOSPHATASE (ACP) Kinetic-spectrophotometric SFBC	40
ALT UV - DAC. Lq ALANINE AMINOTRANSFERASE (ALT/GPT) Kinetic-spectrophotometric IFCC Liquid reagents	200
	1000
ALT-DAC ALANINE AMINOTRANSFERASE (ALT/GPT) Spectrophotometric, end point RAITMAN-FRENKEL METHOD	1200
Albu-DAC. Lq ALBUMIN Spectrophotometric BROMCRESOL GREEN Standard included Ready to use liquid reagent	250
	1000
ALP DEA - DAC ALKALINE PHOSPHATASE (ALP) Kinetic-spectrophotometric DEA BUFFER	100
Amylo D- DAC. Lq α-AMYLASE Kinetic-spectrophotometric DIRECT SUBSTRATE Liquid ready to use reagents	5 x 5 mL
	10 x 15 mL
Amylo-DAC α-AMYLASE Spectrophotometric, end point AMYLOCLASTIC	200
AST UV - DAC. Lq ASPARTATE AMINOTRANSFERASE (AST/GOT) Kinetic-spectrophotometric IFCC Liquid reagents	200
	1000

AST-DAC ASPARTATE AMINOTRANSFERASE (AST/GOT) Spectrophotometric, end point RAITMAN-FRENKEL METHOD	1200
Bili D - DAC. Lq BILIRUBIN (DIRECT) Spectrophotometric DIAZOTIZED SULFANILIC Liquid reagents	100
	500
Bili T - DAC. Lq BILIRUBIN (TOTAL) Spectrophotometric DIAZOTIZED SULFANILIC Liquid reagents	100
	500
Bili DT - DAC. Lq BILIRUBIN (TOTAL AND DIRECT) Spectrophotometric DIAZOTIZED SULFANILIC Liquid reagents	100 + 100
	500 + 500
Bili-DAC-CF BILIRUBIN (TOTAL AND DIRECT) Coffeine method SULFANILIC	DIAZOTIZED 345
Calcium Meth - DAC. Lq CALCIUM Spectrophotometric METHYLTHYMOL BLUE Standard included Liquid reagents	100
	250
Calcium Ars - DAC. Lq CALCIUM Spectrophotometric ARSENAZO III Standard included Ready to use liquid reagent	200
	500
Chol - DAC. Lq CHOLESTEROL Enzymatic-spectrophotometric CHOLESTEROL OXIDASE/PEROXIDASE Standard included Ready to use liquid reagent	200
	500
	1000
Chol HDL - DAC. Lq CHOLESTEROL HDL Precipitation/Enzymatic-spectrophotometric PHOSPHOTUNGSTATE/Mg ⁺² -- CHOLESTEROL OXIDASE/ PEROXIDASE Standard included Ready to use liquid reagent	200

<p>Chol HDL Direct. Lq CHOLESTEROL HDL DIRECT Direct POLYMER/DETERGENT Standard included Liquid reagents</p>	80
<p>Chol HDL - PR CHOLESTEROL HDL PRECIPITATING REAGENT PHOSPHOTUNGSTATE/Mg⁺² Standard included Ready to use liquid reagent</p>	50
<p>Chol LDL Direct. Lq CHOLESTEROL LDL DIRECT Direct POLYMER/DETERGENT Standard included Liquid reagents</p>	80
<p>Chol LDL - PR CHOLESTEROL LDL PRECIPITATING REAGENT POLYVINYL SULPHATE/ POLYETHYLENEGLYCOL Ready to use liquid reagent</p>	20
<p>Cholin - DAC CHOLINESTERASE Kinetic-spectrophotometric BENZOYCHOLINE</p>	20 x 3
	10 x 15
<p>Copper-DAC. Lq Copper Spectrophotometric Dibrom PAESA Standard included Liquid reagents</p>	50
<p>CK - DAC CREATINE KINASE (CK) Kinetic-spectrophotometric IFCC</p>	20 x 2,5
	10 x 15
<p>CK-MB - DAC CREATINE KINASE-MB (CK-MB) Kinetic-spectrophotometric IMMUNOINHIBITION</p>	20 x 2,5
	10 x 10
<p>Crea Kinetic - DAC. Lq CREATININE Kinetic-spectrophotometric ALKALINE PICRATE Standard included Liquid reagents</p>	100
	500
	1000

Crea-DAC CREATININE Spectrophotometric JAFFE Standard included Liquid reagents	260
FRUCTOSAMINE Kinetic-spectrophotometric JOHNSON-BAKER Standard included Ready to use liquid reagent	2 x 50
GGT - DAC. Lq GAMMA-GLUTAMYLTRANSFERASE (g-GT) Kinetic-spectrophotometric IFCC Liquid reagents	50
	200
	500
Gluco - DAC. Lq GLUCOSE Enzymatic-spectrophotometric GLUCOSE OXIDASE/PEROXIDASE Standard included Ready to use liquid reagent	200
	500
	1000
Gluco-DAC GLUCOSE Enzymatic-spectrophotometric GLUCOSE OXIDASE/PEROXIDASE Standard included	200
	1000
Hemo - DAC. Lq HEMOGLOBIN Spectrophotometric ICSH Standard included Liquid reagent, concentrate	1000
Hb Chr - DAC. Lq HEMOGLOBIN CYANIDE-FREE HEMICHROMIC METHOD	3000
Iron Ferr - DAC. Lq IRON Spectrophotometric FERROZINE Standard included Liquid reagents	50
Iron Ch - DAC. Lq IRON Spectrophotometric CHROMAZUROL B Standard included Ready to use liquid reagent	50

TIBC - DAC IRON BINDING CAPACITY AUXILIARY REAGENT	140
LDH - DAC. Lq LACTATE DEHYDROGENASE (LD/LDH) Kinetic-spectrophotometric PYRUVATE Liquid Reagents	50
	200
LIPASE-DAC Kinetic-spectrophotometric Standard included	30
Magne - DAC. Lq MAGNESIUM Spectrophotometric CALMAGITE Standard included Ready to use liquid reagent	50
Phosph - DAC. Lq PHOSPHORUS Spectrophotometric UV PHOSPHOMOLYBDATE/UV Standard included Liquid reagents	170
Potassium-DAC.Lq Potassium Turbidimetry without deproteinization Ready to use liquid reagent	50
Prote - DAC. Lq PROTEIN TOTAL Spectrophotometric BIURET Standard included Ready to use liquid reagent	250
	1000
Prote-DAC PROTEIN TOTAL Spectrophotometric BIURET Standard included	250
	1000
Prot Ure - DAC. Lq PROTEIN (URINE) Spectrophotometric PIRROGALOL RED Standard included Ready to use liquid reagent	50
	250

Sodium-DAC.Lq Sodium Spectrophotometric Standard included	50
Sulfo-DAC PROTEIN (URINE) Spectrophotometric SYLFOSALYCILIC ACID Standard included	200
TG - DAC. Lq TRIGLYCERIDES Enzymatic -spectrophotometric GLYCEROL PHOSPHATE OXIDASE/PEROXIDASE Standard included Ready to use liquid reagent	50
	250
	1000
Thymol-DAC TEST WITH THYMOL Spectrophotometric	500
Ure - DAC. Lq UREA/BUN Enzymatic -spectrophotometric COLOR Standard included Liquid reagents	100
	500
Ure UV - DAC. Lq UREA/BUN Kinetic -spectrophotometric ULTRAVIOLET Standard included Liquid reagents	50
	250
	1000
Uric Acid - DAC. Lq URIC ACID Enzymatic -spectrophotometric URICASE/PEROXIDASE Standard included Liquid reagents	200
	500
	1000
Zinc-DAC. Lq Zinc Spectrophotometric 5-Br-PAPS Standard included Liquid reagents	50

STANDARDS and CALIBRATORS

Product/Description	Total reagent Vol., mL
Gluco-cont GLUCOSE	10
Thymol-DAC-Standard STANDARD FOR TEST WITH THYMOL	125
Ure-cont UREA (THREE CONCENTRATIONS)	3 x 10
Sera L I - DAC ASSAYED CONTROL SERUM LEVEL I	5
Sera L II - DAC ASSAYED CONTROL SERUM LEVEL II	5
Multi St - DAC CALIBRATOR SERUM	5
CK-MB CONTROL SERUM	0,5
Bili St - DAC BILIRUBIN STANDARD	5

DIAGNOSTIC KITS FOR SEROLOGY AND CLINICAL IMMUNOLOGY

Product/Description/Method	Tests Q-ty
ASLO-DAC ANTISTREPTOLYSIN-O LATEX SLIDE TEST	50
	100
CRP-Latex C- REACTIVE PROTEIN LATEX SLIDE TEST	50
	100
IM-DAC INFECTIOUS MONONUCLEOSIS (normal antibodies) HEMAGGLUTINATION SLIDE TEST	50
	100
Hepa Tex-DAC HBsAg LATEX SLIDE TEST	50
	100
RefaTex-DAC RHEUMATOID FACTOR LATEX SLIDE TEST	50
	100
RefaTex-conf RHEUMATOID FACTOR + CONFIRMATION TEST LATEX SLIDE TEST	50 + 50
	100 + 100
SLE-DAC SYSTEMIC LUPUS ERYTHEMATOSUS (DNA antibodies) LATEX SLIDE TEST	50
	100
Waler-Rose - DAC RHEUMATOID FACTOR WAALER-ROSE REACTION HEMAGGLUTINATION SLIDE TEST	50
	100
TPHA SYPHILIS HEMAGGLUTINATION TEST	100
	200
RPR-CARBON SYPHILIS	100

HEMAGGLUTINATION SLIDE TEST	500
PregnaTex-DAC PREGNANCY (hCG) LATEX SLIDE TEST	50
	100

DIAGNOSTIC KITS FOR DETERMINATION OF HEMOSTASIS PARAMETERS

Product/Description	Tests Q-ty
<p>Antithrombin III-DAC</p> <p>ANTITHROMBIN III</p>	60
<p>APTT-DAC</p> <p>ACTIVATED PARTIAL THROMBOPLASTIN TIME</p>	120
<p>Coagulo-cont</p> <p>LYOPHILIZED CONTROL HUMAN PLASMA</p>	3 x 1 mL
<p>Haemostaz-DAC-I</p> <p>ACTIVATED TIME OF RECALCIFICATION</p> <p>ACTIVATED PARTIAL THROMBOPLASTIN TIME (APTT)</p> <p>THROMBIN TIME/PROTHROMBIN TIME WITH PROTAMIN SULFAT</p> <p>PROTHROMBIN TIME, PROTHROMBIN RATIO, INR</p> <p>FIBRINOGEN</p>	60
	60
	100
	60
	60
<p>Haemostaz-DAC-II</p> <p>EUGLOBULIN LYSIS TIME</p> <p>HAGEMAN-DEPENDENT FIBRINOLYSIS (FACTOR XII)</p>	50
	50
<p>Protrombin-DAC-SI</p> <p>PROTHROMBIN TIME</p> <p>PROTHROMBIN INDEX</p> <p>PROTHROMBIN RATIO (PR)</p> <p>INTERNATIONAL NORMALIZED RATIO (INR)</p>	60
	120
<p>Trombin-DAC</p> <p>THROMBIN TIME</p>	100

CLINICAL IMMUNOLOGY

	Total reagent Vol., mL
CIC-DAC Circulating immune complexes Spectrophotometric Polyethyleneglycol	435
NBT-DAC Bactericidal activity of phagocytes Spectrophotometric Nitroterazolium blue reduction test	304

BLOOD GROUPING REAGENTS

Product	Total reagent Vol., mL
Anti-A	10
Anti-B	10
Anti-AB	10
Anti-D-IgM	10

RAPID TESTS "ONE STEP"

Product/Description	Tests Q-ty
Pregnancy test (hCG), urine	1
Ovulation test (LH), urine	1
HBsAg test (HBs antigen), serum/plasma	1
HBsAg test (HBs antigen), whole blood	1
Anti-HBs test (Hepatitis B Surface antibodies), serum/plasma	1
Anti-HBs test (Hepatitis B Surface antibodies), whole blood	1
HBeAg/HBsAg test, serum/plasma	1
Anti-HCV test (Hepatitis C Virus antibodies), serum/whole blood	1
Anti-HIV 1/2 test, serum/whole blood	1
Anti-HIV 1/2 test, serum/whole blood	25
Anti-H.Pylori test, serum/whole blood	1
Anti-Treponema Pallidum test (syphilis), serum/whole blood	1
Anti-Treponema Pallidum test (syphilis), serum/whole blood	25
Chlamidia Antigen test, serum/whole blood	1
AFP test (alpha-fetoprotein), serum/plasma	1
AFP test (alpha-fetoprotein), serum/plasma	25
CEA test (carcinoembryonic antigen), serum/plasma	1
PSA test (prostate specific antigen), whole blood	1
PSA test (prostate specific antigen), serum/plasma	1
CEA/AFP test, serum/plasma	1
MOP test (Opiates, Morphine, Heroin), urine	1
MET test (Methamphetamine), urine	1
AMP test (Amphetamine), urine	1
COC test (Cocaine), urine	1
THC test (Marijuana), urine	1
Alcodiagnostic (alcohol detection in saliva)	25

Urine strips

Product		Tests quantity
DAC-1G (Glucose)		100
DAC-1K (Ketones)		100
DAC-1P (Protein)		100
DAC-1B (Blood)		100
DAC-2GK (Glucose, Ketones)		100
DAC-2GP (Glucose, Protein)		100
DAC-2BU (Bilirubin, Urobilinogen)		100
DAC-3 (Glucose, Protein, pH)		100
DAC-4B (Glucose, Protein, pH, Blood)		100
DAC-4S (Glucose, Protein, pH, Specific Gravity)		100
DAC-5 (Glucose, Protein, pH, Blood, Ketones)		100
DAC-7L (Glucose, Protein, pH, Blood, Nitrite, Specific Gravity, Leucocytes)		100
DAC-8N (Glucose, Protein, pH, Blood, Ketones, Bilirubin, Urobilinogen, Nitrite)		100
DAC-8S (Glucose, Protein, pH, Blood, Ketones, Bilirubin, Urobilinogen, Specific Gravity)		100
DAC-9 (Glucose, Protein, pH, Blood, Ketones, Bilirubin, Urobilinogen, Nitrite, Specific Gravity)		100
DAC-10 (Glucose, Protein, pH, Blood, Ketones, Bilirubin, Urobilinogen, Nitrite, Specific Gravity, Leucocytes)		100
DAC-11 (Glucose, Protein, pH, Blood, Ketones, Bilirubin, Urobilinogen, Nitrite, Specific Gravity, Leucocytes, Ascorbic Acid)		100